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Cytes Biotechnologies供給

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<https://www.cytesbiotechnologies.com/>

Cryopreserved Human Hepatocytes Available

Use the **filters** below to find the best lot for your needs.
Click on the "more info" column to get the **characterization datasheet**.
Use the **export option** on the right hand side of the filters to export the data
HuCPM: Certified **M**etabolism studies / **HuCPI**: Certified **I**nduction studies / **HuCPT**: Certified **T**ransporter studies

Plateable Non Plateable

Metabolism Certified: Induction Certified: Transporter Certified: Stock: 330 Cells/vial post thaw: 20 Viability post thawing: 100%

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	Product Referer	More Info ↑	Lot ↑	Stock	Viability	Cells per vial (M)	Age	BMI	Gender	CYP1A n-fold	CYP2B n-fold	CYP3A n-fold
1	HuCPM+	I-M	https://goo.gl/zf BHum16061	28	87.1%	7.9	74	27.34	Male	21.79	2.39	2.66
2			https://goo.gl/8d6 BHU416043	112	92.7%	6.7	74	26.9	Female	26.88	6.9	11.76

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<付属データ>

- Donor Demographics
- Post Thaw information
- PHASE I: P450 activities expressed

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- PHASE I: P450 activity & mRNA fold-induction
- PHASE II: UGTs & SULT activities expressed

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87%*

*Post Thaw information
2017年調査

Certification date: 1st February 2019

ANALYSIS CERTIFICATE

Lot#: BHum16016

DONOR DEMOGRAPHICS

Species	Sex	Race	Age	BMI	Smoker	Alcohol Use	Drug Use
Human	Male	Caucasian	70	26	No	No	No

Pathology	Medications	Serological Data
Colorectal cancer hepatic metastases	Paroxetine, duodart	Tested negative less than 3 months before surgery

Patient informed consent was obtained. The donor was serologically tested negative for following infectious diseases: HIV, Hepatitis B and C. Donor medical history were also examined prior to accepting this donor.

For in vitro use only, not to be used for clinical application. Products distributed by Cytes Biotechnologies may contain human material that should be treated as potentially hazardous.

CHARACTERIZATION FOR PLATEABLE CELLS

Lot#: BHum16016

Post Thaw Lot information

Number of viable cells/vial:

9.5 ± 1.3 × 10⁶ (n=2)

Post-thaw viability (%):

90.3 ± 2.7 (n=2)

Monolayer assessment*



Plateability: YES

Seeding density in 24 well recommended:

0.3 × 10⁶ cells/well in 0.5 ml

Cell morphology 24h

Cell morphology 96h

• PHASE I: CYP ACTIVITIES EXPRESSED IN pmol/min/mg protein (mean ± SD)

Enzyme	Induction (Specific Activity)			
	Basal Activity (pmol/min/mg protein) on day 1	Basal Activity (pmol/min/mg protein) on day 4	Induced Activity (pmol/min/mg protein) on day 4	n-Fold induction
CYP1A2	11.52 ± 1.69	1.95 ± 0.06	341.66 ± 84.99	175.08
CYP2B6	13.19 ± 0.31	8.09 ± 0.81	74.83 ± 5.34	9.25
CYP3A4	13.73 ± 0.54	0.94 ± 0.12	20.71 ± 5.20	21.93

Cryopreserved human hepatocytes were thawed and plated in 24well collagen I coated plates. Cells were overlaid (Corning) in Human Hepatocyte Maintenance Medium at first medium change at day of thawing. Treatment (n=2 per vehicle control [0.15% (v/v) DMSO] or inducers (Rifampicin, β-Naphthoflavone and Phenobarbital) began 1-day post-plating and continued for 72 hours. At the end of induction, monolayers were rinsed with PBS and incubated with probe substrate culture media. See Table 1 for information on each probe substrate. Metabolites were quantified by LC-MS and normalized to the vehicle control. The fold induction was calculated by dividing the induced activity by the vehicle basal activity on the same day.

• PHASE I: CYP450 mRNA induction

CYP (mRNA)	n-Fold Induction
CYP1A2	11
CYP2B6	22
CYP3A4	191

Cryopreserved human hepatocytes were thawed, plated in 24well collagen I coated plates in Hepatocyte Plating Medium overlaid with Matrigel® (Corning) in Human Hepatocyte Maintenance Medium at first medium change at day of thawing. Treatment (n=2 per compound) with vehicle control [0.15% (v/v) DMSO] or inducers (Rifampicin, β-Naphthoflavone and Phenobarbital) began 1-day post-plating and continued for 72 hours. At the end of the period, RNA was isolated for mRNA analysis.

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